

## **Pulse Oximetry during Low Perfusion Caused By Emerging Pneumonia and Sepsis in Rabbits.**

Hummler H.D., Pohlandt F., Franz A.R. *Crit Care Med.* 2002 Nov;30(11):2501-8.

### **Objective**

This study tested the effects of low perfusion caused by emerging sepsis on the reliability of a new pulse oximetry technology (Masimo SET; IVY 405T) compared with a standard pulse oximeter (Nellcor N-200).

### **Methods**

*Design:* Randomized trial. *Setting:* University animal research facility. *Subjects:* Twenty-six anesthetized, ventilated (FiO<sub>2</sub>, 1.0), adult rabbits. *Interventions:* Pneumonia/sepsis was induced by tracheal instillation of Escherichia coli. Oxygen saturation was measured by pulse oximetry (SpO<sub>2</sub>) and recorded continuously until death. Arterial oxygen saturation (Sao<sub>2</sub>) was measured hourly by oximetry and whenever SpO<sub>2</sub> dropped to  $\leq 95\%$ , or whenever a difference of  $\geq 5\%$  between devices occurred. SpO<sub>2</sub> sensors were positioned at both forelegs and switched hourly.

### **Results**

The total time of signal loss was longer with the N-200 vs. the IVY: 65 (4-299) min vs. 7 (0-97) min [median (range)],  $p < 0.001$ . Signal loss was more prevalent during the first 80% of the experimental time with the N-200 compared with the IVY. Nineteen of 26 animals had a total of 62 episodes of a falsely low SpO<sub>2</sub> value with either one of the two devices associated with hemodynamic deterioration. Median bias (SpO<sub>2</sub> - SaO<sub>2</sub>) was small, but variability of bias values increased toward the end of the experimental time with both devices.

### **Conclusions**

The pulse oximeter equipped with Masimo SET was less prone to signal loss than the standard pulse oximeter in this sepsis model. Episodes of falsely low SpO<sub>2</sub> readings may occur, and deviation of SpO<sub>2</sub> from SaO<sub>2</sub> may be increased with deteriorating hemodynamics with both devices.